

DigiTrak XT Recorder— smart and innovative

On-screen lead map and EASI lead placement for artifact-free

12-lead standard on all recorders

Pacemaker detection with sensitivity adjustment

Lightest recorder with large onboard display, just 2.2 oz (62 g)

Extra-slim design for patient comfort

Swivel belt clip and pouch for easy wearability

Recording up to 7 days on a single AAA battery

Water-resistant case and 2-year warranty

Non-removable memory prevents accidental data loss and enhances data security




Automatic Power On self-test with onscreen status indicator guarantees reliable recording

Simple navigation buttons with tactile response for fast setup

Large, raised Event button easy for patients to find by feel

This unit is configured for:
 First Name: JOE
 Last Name: SMITH
 Patient ID: 0123456789
 Press < and ← to accept

Streamline your workflow by pre-configuring and previewing patient demographics.



Large display for each channel ensures artifact-free



Fast downloads typically take less than 90 seconds. Our dual dock accepts DigiTrak XT and DigiTrak Plus recorders.