

Uncompromising Clinical and Technical Holter Performance.

The thin and lightweight Chroma Holter recorder provides superior pacemaker detection with adjustable sensitivity for the latest pacemaker technology. Add advanced analysis capabilities with easy to use, feature-rich HolterCare™ software.

FEATURES

- **Records up to 30 days of continuous, 3-channel, ECG data**
- **Data Transfer Options:** Choice of direct USB data download or removable data memory card
- **Cable Detection:** Smart cable detection between 5 and 7 lead, with display of the appropriate placement based on the cable attached
- **Direct Demographics Downloads:** Downloads patient demographics into Chroma via automatic recorder setup through direct HL7 interface
- **Adjustable Data Sample Rates:** Adjustable data sample rates from 128 to 1024 samples per second with 8, 10 and 12 bit resolution for the most demanding clinical requirements and analysis
- **Notifications:** Audible and visual patient lead disconnect and medication/activity reminders
- **Hook-up Guidance:** On-screen hook-up chart with visual guidance for proper electrode placement and lead disconnect indication
- **Multiple user-selectable languages:** Includes English, Spanish, French, Portuguese, Italian, German, Russian, Dutch, Finnish and Danish
- **Customizable:** Fully customizable recorder interface and controls
- **Convenient:** Includes neoprene pouch with belt clip



BENEFITS

- Preserves data integrity through superior ECG baseline correction for cleaner data
- Downloads ECG data in less than 15 seconds, with automated analysis in 60 seconds
- Eliminates paper diary transcriptions through easy-to-use, color coded, digital diary
- Easier to use thanks to a large, bright, color OLED display
- Confirms signal quality through visual ECG display

Chroma: Options for Every Need

A variety of optional Chroma features and additional products will further enable you to maximize diagnostic and operational efficiencies.

MORE ADVANCED DATA FOR COMPLETE INSIGHT

- Multi-day Holter analysis
- Simple and reliable pacemaker analysis
- Statistical and spectral heart rate variability through HolterCare

QUICK DATA ACCESS AND TRANSFER

- EZLink allows automated transfer and analysis via the Internet

ENSURED COMPLIANCE

- HIPAA-compliant software with HL7 Interface supports EHR connectivity

SINGLE PLATFORM SUITE ACCESS

- Chroma and HolterCare are part of the CardioView™ Dx Suite, the industry's most comprehensive single-platform, fully integrated product suite for IDTFs, hospitals and physicians
- Uniquely facilitates the management of multiple ambulatory ECG services and remote cardiac rhythm device diagnostics while reducing costs and maximizing efficiencies


CardioView Dx™ DIAGNOSTIC SUITE

Chroma & HolterCare are part of the CardioView™ Dx Suite
LIVE AMBULATORY TELEMETRY • ONEVIEW CRM • HOLTER • ABP

The HolterCare™ Advantage

MAINTAIN COMPLETE CONTROL OVER YOUR DATA

- One mouse click provides vital information, including all arrhythmic and ischemic episodes, which can be visually confirmed and edited on a beat-by-beat basis
- HolterCare provides comprehensive arrhythmia editing techniques, including multi-level template editing, along with superimposition and page mode review
- Preview data set and set the correct channels of analysis before editing, allowing analysis of the best channels, thus yielding the best results most efficiently



ROBUST DATA WITH MEANINGFUL PRESENTATION

- Color-differentiated beats, bi-directional split-screen superimposition, page mode and colored graphic trends help the user easily locate areas of diagnostic interest
- FDA XML reporting capabilities present only the data the FDA needs to see, greatly improving efficiencies
- Automatic a-fib identification and burden calculation utilizing industry recognized algorithm